Profi

Original

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Operating Instructions

Mode d'emploi

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D
E
F
I
ES

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Original

Operating Instructions

We are pleased that you decided to buy a highly developed piece of equipment from SCHICK and would like to wish you every success when working with your new  **Profi** .

We wrote these operating instructions to enable you to get accustomed to your new piece of equipment and to provide you with the correct operating and maintenance instructions.

Index	Page
1. Safety Instructions.....	3
Caution!	
2. Range of applications.....	3
3. Summary List of content.....	4
4. Equipment/Optionals.....	5
5. Initial start up and installation.....	6
Switching on - and off.....	7
6. Tool change on motorhandpiece.....	8
7. Operation.....	8
8. Cleaning and maintenance/.....	
Removing the chuck.....	9
9. Errors.....	10
10. Technical data.....	10
11. Declaration of conformity.....	11

1. Safety information

- 1.1 Ascertain that your mains supply coincides with the data on the rating plate.
- 1.2  **Profi**- units are not suitable for the following applications:
- in areas where there is a risk of explosion
 - for medical applications
 - for working on moist materials
- 1.3 Ensure that all regulatory requirements are observed during use
- always wear protective glasses
 - provide enough light at the working place
 - use dust suction
- 1.4 Under no circumstances should the motor handpiece be cleaned with compressed air.
- 1.5 Before putting the handpiece down, always insert a rotary instrument or the pin, supplied with the unit, into chuck.

Caution!

- When using rotary instruments, do not exceed the maximum speeds laid down by their manufacturer.
- When operation with left-hand rotation, the collet may loosen when there is a high load.
- Repairs and other technical procedures must only be carried out by suitable qualified personnel, authorized by SCHICK.
- SCHICK do not guarantee the  **Profi** unit should it not have been used in accordance with the operating instructions.

2. Range of applications

The  **Profi** unit is designed for universal use in dental laboratories when trimming crowns and bridges, respectively acrylic and light chrome cobalt dentures. With its speed range from 1.000 - 40.000 rpm the SCHICK  **Profi** - unit allows to work all dental materials.

The motorhandpiece is ergonomical, small shaped and runs smoothly, thus eliminating user-fatigue and making it economical to use over long periods. The progressive processor controller constantly monitors the unit for overloading. A second thermal fuse, which is independent of the computer, provides further protection. Thus the  **Profi** - unit offers a maximum of safety.

Conditions of environment:

- interior 5° - 40° C.
- up to 2,000 meter over sea level

Categorie of overvoltage: II

Grade of pollution: 2

WEEE-Reg.-Nr. DE 78620387

3. Summary - List of contents

(1)


(2)

(3)

(4)

(5a)

(5b)

(6)

(7) (8)

(9)

kneeverision:

Profi SK complete		Art.Nr.	9340
Profi knee controller	(1)		9345
motorhandpiece with cable	(3)		9310
handpiecerack	(4)		9127
with two keys for changing the chuck			
AC adaptor	(5a)		9102
velcro	(5b)		9130
mains cable	(6)		2160
suspension strip	(7)		9103
screws - two pieces-	(8)		3170
signaltransmitter for suction	(9)		9060

foot control:

Profi SF complete		Art.Nr.:	9350
Profi foot controller	(2)		9355
motorhandpiece with cable	(3)		9310
handpiecerack	(4)		9127
with two keys for changing the chuck			
AC adaptor	(5a)		9102
velcro	(5b)		9130
mains cable	(6)		2160
signaltransmitter for suction	(9)		9060

4. Accessories/optionals

chuck Ø 2,35mm	6642
chuck key	4115
counterstay wrench	6223
motor cable complete	9225

riveting hammer
Art.no. 1850

twin riveting hammer
Art.no. 1860


riveting hammer with chisel
Art.no. 1850/1


5. Initial start-up and installation

5.1 Installation of the suspension strip for knee control unit or work bench control unit

1.

2.


- 1 The suspension strip (7) is supplied with a template which acts as an aid for correct positioning of the **Profi** knee control unit. To do this, the mounting bracket is positioned on the work bench as shown in figure 1, and secured using the screws supplied (8). After installation, the template can simply be detached from the bracket.
2. For installation of the knee control unit, the cavity at the rear of the unit is placed on the bracket and pushed back until it clicks into place (Fig.2).

5.2 Fixing the power supply with the velcro

The delivered velcro can be used to attach the power supply under the workbench (Fig.2)

1.

2.


1. Remove one of the protective foils of the velcro and attach it to the power supply.
2. Remove second protective foil and stick the power supply to a clean and even surface. Please avoid permanent load to the power supply by the cables.

5.3 Initial start-up

Connect the network cable (6) to the power supply (5a). Connect the motor hand piece to the 'motor hand piece 1' socket.

Plug the network adaptor into the power supply socket on the control unit.

To switch on the device, activate the operating component (knee pad, foot lever, etc.) once.


5.4 Switching on and off:

There is no main switch for the **Q-Profi** control unit. If the hand piece is not used for more than 60 minutes, it switches off automatically. By simply touching the control device (knee pad /foot pedal), the device is switched on again.

 Basically, first attach the handpiece before connecting the mains cable and switching the control mechanism on.


5.5 Connection of the signal transmitter for automatic switching on of a workplace suction device

In order to operate the workplace system together with a workplace suction system with automatic function, the signal transmitter for the suction signal must be connected. This is, as shown in the illustration, inserted in the suction signal jack on the control unit. At the opposite end, in the appropriate connector socket on the workplace suction device.

When operation **Q-Profi** together with a KaVo SMARTair suction or a Zubler suction AP 501, the optional cable Art.No.: 9229 can be used as signal connection for controlling the suction.

6. Tool change on motorhandpiece


By turning the gripsleeve on the hand piece, the chuck can be opened or closed.

- ⚠ Only carry out a tool change when the motor is switched off! To optimise the precision and durability of the chuck, the tool or the factory-supplied protective pin must always be in a fixed position, even when not in use..
- ⚠ Always insert tools as far as possible into the chuck in order to achieve maximum holding power!

7. Operation

There is continuously variable speed control of the handpiece up to the maximum preset motor speed using the respective control unit (see diagram below).


The knee controller is equipped with an additional potentiometer for setting a speed limit. Set the potentiometer to the desired maximum speed. By using the knee lever the speed can now be set variously up to the pre-set speed.

- ⚠ Please always observe the maximum permitted speed for your tools!!

8. Cleaning and Maintenance / Removing the chuck

The **QProfi** motor handpiece is designed for maximum durability and therefore, the chuck and the tip should be removed and cleaned from time to time.

8.1 Removal of chuck on the **QProfi** handpiece

1. Remove tool and unscrew tip of handpiece
2. Remove cap from motor and pull out the plug
3. Place the chuck key (Art. no. 4115) on the open chuck, use a number 6223 counterstay wrench to hold the motor end of the spindle
4. Remove the chuck from the shaft

After thoroughly cleaning the chuck, before replacing, it should be slightly greased on the outer surface (e.g., with Vaseline).
Only fasten the chuck lightly when reinstalling! The chuck will fix itself firmly as part of its normal operation.

To optimise the precision and durability of the chuck, a tool or the factory supplied protective pin must always be fixed, even when not in use.

Never clean the handpiece with compressed air!

The required chuck keys are located in the underside of the handpiece holder.

8.2. Removing/changing the motor cable (see picture above)

1. Unscrew the cap (2) from the motor
2. Loosen the cable by pulling out the plug-in seal (1)

9. Errors

If the handpiece is overloaded or blocked, the control unit switches off automatically to prevent further damages. Make sure the tool can rotate freely in the handpiece and after releasing foot lever or knee lever back to neutral position the unit can continue operation.

If an error cannot be removed following the above description, please contact an authorised service partner or Schick directly.

10. Technical data

dimensions handpiece:

length: 158 mm
diameter max.: 30 mm
weight : 240 g

dimensions control unit:

	knee version	foot version
width:	95 mm	155 mm
height:	203 mm	75 mm
depth:	210 mm	220 mm
weight:	902 g	1,82 kg

operating voltage: AC 100 - 240 Volts
Frequency: 50/60 Hz
output: 160 Watt

speed range: 1.000 - 40.000 min
torque: max: 7,5 Ncm
protection class: I

The oscillation total value during operation is below $2,5 \text{ m/s}^2$ (EN 28662)².

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declare herewith, that the products

- Profi** - consisting of
- Profi** - handpiece 9310 in combination with
- Profi** - control unit 9345 or 9355

are in conformity with the following provisions of Directive:

- 2001/95/EG (general product safety)
- 2006/42/EG (machinery directive)
- 2006/95/EG (low voltage directive)
- 2004/108/EG (EMV directive)

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Schemmerhofen, March 2013

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Subject to technical modifications

This unit complies with the current VDE (German association of electrical technicians) regulations concerning safety and suppression.

These operating instructions should be readily accessible and are best kept close to the unit itself.

We would like to take this opportunity to advise you that a proper repair service and suitable qualified personnel are required for such highly developed technical equipment. SCHICK guarantees to carry out perfect repairs using original spare parts.

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Schick GmbH